

himac
EPPENDORF GROUP



eppendorf
welcomes Himac as a
new member of the
Eppendorf Group

CP-NX SERIES

Ultracentrifuge

Ultracentrifuge



himac
EPPENDORF GROUP

CP-NX SERIES

CP100NX / CP90NX / CP80NX

Truly User-friendly Design

LED Indicator – Operation status is at a glance !

LED Indicator of 47cm length is mounted in front of table. This indicates its operating status by selectable light in colors, luminous patterns, and levels of brightness. So operating status is at a glance with this LED Indicator. These are :

Status of operating:

Stop, acceleration, running at set speed, deceleration, vacuum holding, economy mode, acceleration in zonal mode, deceleration in zonal mode, or alert.

Color:

Blue, red, green, light blue, yellow, white or pink



Brightness:

Three levels

Luminous pattern:

Solid lighting, slow and fast blinking, or fading

* Actual color of LED Indicator may be different from the photos on this information, brochure and website due to printing conditions.



Touch-sensitive LCD panel with intelligible screen design

The control panel, 6.5-inch size touch-sensitive color LCD panel, is located in a front part of the machine. You do not need to stretch your arm to operate the centrifuge. The intelligible screen designs, such as a wide indication of basic parameters, many icons for selecting function, folder colors with each symbol for program management, offer user-friendly operation, like operating a smart phone, to all users.



Multilingual Display

You can select a language from 11 languages, English, German, French, Spanish, Italian, Portuguese, Netherlands, Russian, Chinese (simplified), Korean and Japanese.



Low table height with smooth sliding door for easy loading / unloading the rotor

The table height of the CP-NX is 863mm. The low table height enables users to load / unload the rotor easier. The rotor is automatically locked by the centrifugal force once the rotation is started by the "self-locking rotor system". So it is not necessary to fix the rotor by screwing the rotor onto the drive shaft or push the button to lock/unlock the rotor any more. In addition, the improved door structure design based on the latest structure analysis technology makes the door thickness 5mm thinner than the one of our former model. The lighter door offers the same protection as the before and makes door-sliding movement (back and forth) smoother than before.






Advanced Technologies, Features & Functions

Automatic Rotor Life Management System

Automatic Rotor Life Management System, RLM automatically manages an operating record of a rotor to accumulate totals of run times and operating hours by each rotor.

RLM rotor has a magnetic memory on its rotor itself as RLM Adapter. This system is original only for himac; the memory is read out and over-written by RLM sensor with a himac ultracentrifuge.

himac ultracentrifuge reads out a record of run times and operating hours of rotor from its memory during acceleration to accumulate the both times and hours by RLM adapter and sensor. While the rotor decelerates, the totals of the run times and hours are accumulated and updated by centrifuge, the latest accumulated record are overwritten into memory on RLM adapter of the rotor.




Though a RLM rotor installed on and run by himac's another centrifuge with RLM system, an operating record of a RLM rotor is kept as the rotor's latest of the accumulated running run times and hours. Because the accumulated record is always updated and overwritten into memory of RLM rotor.

Longer service of Rotor

Automatic Rotor Life Management System tracks a rotor running precisely in times and hours. It is more accurate than management manually by the rotor log book and its conventional calculation.

When a rotor runs lower speed than the maximum, and furthermore, its running time is less than an hour, this RLM system evaluates its accumulated record to be less than actual running in time and hour, corresponding its smaller load on the rotor. Consequently, the rotor's expected service is much longer in time and hour by this system, compared with management by log book.




If a rotor speed is lowered from 100,000 rpm to 80,000 rpm by 20% less, its expected services are increased from 5,000 runs to 10,000 runs as a double.

If a rotor speed is lowered from 100,000 rpm to 50,000 rpm as a half, its expected service hours are increased from 10,000 hours to 20,000 hours as a double.

Balancing by non-contact imbalance sensor

Easy to balance all the sample by balancing them within 5 mm among the all tubes and bottles to be set in a rotor. Non-contact imbalance sensor always monitors of both the rotor and shaft vibrating. When the both vibrate unusually, the centrifuge automatically stops the rotor.

Remark: This feature does not applied with rotor P21A2.



Economy mode* for energy saving


Stand-by electricity consumption reduced to a half **, compared with our former model, by deeming LCD backlight, stopping cool system and its fans when the control panel untouched over the set period in 1 to 180 minutes.

* This economy mode works while drive unit and vacuum pump are stopped.
** Under in-house test condition



Operating History : 5,120 logs

5,120 run histories recorded automatically and reuseable in user program. A run history exportable through USB interface in CSV format.



Programmed Operation : 1,000 programs

1,000 use-programable operations with file folders in 4 selective color.



User Administration with User Lockout

50 users can be classified into 3 levels in access: Administrator, Supervisor, and User, secured by passwords to enter log-in system for each access level and its functions.



Quiet and Fast

Quiet operation

Operating noise level 51 dB(A) under running at maximum speed. As of Jan. 2014, under in-house test condition

Quick start-up

6 seconds and ready to operate system after power-on.

Faster vacuum

Evacuation time 15 minutes to reach high vacuum. The CP-NX series has reduced 5 minutes to reach the high vacuum, compared with our former model.

As of Jan. 2014, by our investigation
No condensation in a rotor chamber before starting vacuum



Intelligible Screen Design

RUN Screen



Program Button
To set, select and manage programs for programmed operation



Vacuum Button
To start and stop vacuum, also shows vacuum level



RCF Button
To display and set RCF



Start Button
To start operation



ω²T Button
To set ω²T operation



Stop Button
To stop operation



RTC Button
To set timer (RTC : Real-Time Control) operation

MENU Screen



Run History
Automatically record operation log up to 5,120 runs for past runs.



Customize
To enter Customizing Screen



Rotor Catalog
To browse available rotors, their specifications and accessories



Manager (Admin)
To enter Administration Screen



Zonal Setting
To switch between normal and zonal operation



Defrost Function
To defrost rotor chamber



ID/Contact
To key-in instrument's ID and service contact information

CUSTOM Screen



Zoom
To enlarge Speed and Time indication on RUN screen during the running



Schedule
To record usage schedule up to 40 schedules



Stop Signal
To select sound of stop signal



Date/Time
To set date and time



Volume
To adjust sound volume of the stop signal



LED Indicator
To set LED Indicator for each operating status (colors, luminous patterns, brightness)



Backlight
To adjust backlight of LCD panel



Economy Mode
To set economy mode for energy saving



Language
To select display language from 11 languages (English, German, French, Spanish, Italian, Portuguese, Nederlands, Russian, Chinese, Korean and Japanese)

ADMIN Screen



User Management
To register, manage and delete users



Actual Run Timer
To set actual run timer (excluding acceleration time from run time)



User Lockout
To limit access of users



Vacuum Level
To set vacuum level of acceleration from vacuum stand-by to set speed



Rotor Management
To manage total number of run and run time of registered rotors



Zonal Speed
To set zonal speed between 2,000 to 3,000rpm with 100rpm increments



Run Time Indication
To change run time indication between Elapsed Time or Remaining Time



LAN Communications
To connect or disconnect LAN communications



Data Communication and Log Management Software

USB (device) x 1, USB (host) x 1 and LAN x 1 are equipped as an interface for the data communication. Operating histories, up to 5,120 logs, can be exported in CSV format through USB (host) port.



himac LogManager ver. 5.0 for windows® (network edition) – optional log management software


himac LogManager ver.5.0 for Windows® (network edition) is a useful and convenient software to manage real-time operating log of himac CP-NX series, CP-WX series, CS150NX, CS150FNX, CR30NX, CR22N and CR7N. Maximum 16 units of the above-mentioned centrifuges can be registered to the software and monitored by the software at the same time. The software is installed into the PC and data communication between the centrifuges and the PC is done via LAN or Ethernet. You can easily establish the network configuration and relocation with commercially available LAN devices. (optional LAN board is required for himac CP-WX series, CS150NX, CS150FNX, CR22N and CR7N.) It means you can manage the operation log at the different location from the installation site of the centrifuges.

The software is real-time log management software; recording interval period is selectable from 10 seconds to 5 minutes. Of course, the software supports U.S. FDA 21 CFR Part 11, following functions are available;

- Digital Signature
- Audit Trail
- Encrypted Data Files


Required Operating System :
Windows®8 and Windows®10

Windows®, Windows®8 and Windows®10 are registered trademarks of Microsoft Corporation in the United States and other countries. Ethernet is a registered trademark of Xerox Corporation in the United States and other countries.



Remote monitoring System himac View

Remote monitoring application "himac View" is available for CP-NX series, CS150NX, CS150FNX and CR30NX. "himac View" monitors, operates or stops a centrifuge with iPhone™ or Android™ Phone via Wi-Fi®. Download "himac View" from Apple™ store or Google play™. "himac View" monitors up to 16 units of CP-NX series, CS150NX, CS150FNX or CR30NX in the LAN access network.



Alert screen



Running screen



Acceleration screen

*Apple and iPhone are trademarks of Apple Inc.
Google, Android, and Google Play are trademarks of Google Inc.
Wi-Fi® is a registered trademark of Wi-Fi Alliance.
Windows® is a registered trademark of Microsoft Corporation in the United States and/or other countries.




Simulation Software "himac ASSIST" – centrifugation simulation software

himac ASSIST is simulation and calculation software of centrifugal condition and can be installed in your Window®-based PC. CD-ROM of the himac ASSIST is included in the standard accessories. You can simulate whether the centrifugal condition is appropriate before the centrifugation, also can simulate optimal centrifugal condition of a sample, based on your CP-NX series ultracentrifuge and rotor.

himac ASSIST has following functions;

- Calculations of K-factor and pelleting time
- Calculation of the allowable rpm with high-density liquid
- Rate zonal centrifugation simulation
- Isopycnic centrifugation simulation
- Solvent concentration conversion
- Mutual conversion of molecular parameters
- Rotor database

Example : Simulation of isopycnic separation




Plasmid DNA separation
with P100VT vertical rotor
CsCl+EtBr : initial density = 1.55g/ml
Speed : 100,000rpm



Biosafety Option

The biosafety is always priority in laboratories. In order to prevent from exhausting bio-hazardous sample into a room, a Micro filter can be assembled in a vacuum line at option. For details, please contact us.





New Rotor

■ Swing Bucket Rotors



Top loading of Titanium rotor



Bucket made of Titanium

Top Loading

Max Speed 32,000 rpm

Application :

- Separation of virus particle
- Preparation of exosome analysis
- Density gradient centrifugation of intermediate volume samples

P32ST



Specifications

| Model | P32ST | P32ST2 |
|------------------------|---|---|
| Max. Speed | 32,000 rpm | 32,000 rpm |
| Max. RCF | 180,000 xg | 193,000 xg |
| K-factor at max speed | 198 | 216 |
| Nominal Capacity | 40 mL x 6 tubes | 16 mL x 6 tubes |
| Loading | Top loading | Top loading |
| Rotor Weight | 7.1 kg | 7.2 kg |
| Accessory included | 100 tubes of 40PA Tube 6pairs of titanium buckets and caps | 100 tubes of 16PA Tube 6pairs of titanium buckets and caps |
| Material No. | 5720214003 | 5720214004 |
| Applicable centrifuges | CP-NX, CP-WX and CP-MX series* | CP-NX, CP-WX and CP-MX series* |

P32ST2



■ Fixed Angled Rotor

P21A2



Application for large volume centrifugation
Nominal capacity 230 mL x 6 bottles



Specifications

| | |
|------------------------|---------------------------------|
| Model | P21A2 |
| Max. Speed | 21,000 rpm |
| Max. RCF | 71,000 xg |
| K-factor at max speed | 486 |
| Nominal Capacity | 230mL x 6tubes |
| Angle of tube cavity | 26 degree |
| Rotor Weight | 11.4 kg |
| Accessory included | 6 pairs of 230PA Bottle and Cap |
| Material No. | 5720211205 |
| Applicable centrifuges | CP-NX, CP-WX and CP-MX series* |

Note: Visual balancing of sample surface difference is within 2mm between the all tubes for rotor P21A2.

Exclusive 230PA Bottle for P21A2



- Thick walled for high centrifugal force
- Made of PPCO, Polypropylene Copolymer for high resistance against cemicals
 - Inner cap made of POM, Polyacetal

| Description | Quantity | Material No. | Remarks |
|--------------|------------|--------------|---|
| 230PA Bottle | 10 bottles | 5720411169 | Excluding cap. |
| AL Cap (2) | 2 caps | 5720411170 | Exclusive cap for rotor P21A2 and 230PA Bottle. Including O-ring. |
| O-Ring | 10 rings | 5720411174 | As a spare of AL Cap (2) |

Remarks :

* To use these new rotors with himac former ultracentrifuges, CP-WX or CP-MX series, ROM in CP-WX or CP-MX series may be exchanged based on its manufacturing number. Please inform us model name and the number of your ultracentrifuge to check if the ROM should be exchanged or not before placing the rotor.




Rotor List

| Type | Model | Material No. (w/RLM Adapter) | Max. Speed (rpm) | Max. RCF (xg) | Nominal Capacity (ml x tubes) | K-factor | Main Purpose |
|-----------------------------|------------|---------------------------------|---------------------|------------------|-------------------------------------|----------|--|
| Fixed Angle Rotor | P100AT2* | 5720211101 | 100,000 | 803,000 | 6.5 x 8 | 18 | Separation of microscopic particles and lipoprotein |
| | P90AT* | 5720211102 | 90,000 | 700,000 | 12 x 8 | 25 | Separation of cell organelles, Plasmid DNA and RNA Rapid separation of intermediate volume sample |
| | P80AT* | 5720211103 | 80,000 | 615,000 | 12 x 8 | 27 | |
| | P70AT2* | 5720211104 | 70,000 | 452,000 | 12 x 12 | 36 | |
| | P70AT* | 5720211105 | 70,000 | 505,000 | 40 x 8 | 44 | Rapid separation of intermediate-to-large volume sample |
| | P65A*† | 5720211201 | 65,000 | 370,000 | 12 x 10 | 48 | Separation of cell organelles, Plasmid DNA and RNA |
| | P50AT4 † | 5720211202 | 50,000 | 316,000** | 6.5 x 44 | 32** | Separation of lipoprotein |
| | P50AT2* | 5720211107 | 50,000 | 303,000 | 40 x 12 | 70 | Rapid separation of intermediate-to-large volume sample |
| | P50A3 | 5720211108 | 50,000 | 252,000 | 1.5 x 24 | 33 | Rapid pelleting micro volume samples |
| | P45AT* | 5720211109 | 45,000 | 235,000 | 94 x 6 | 130 | Rapid separation of large volume samples |
| | P42AT † | 5720211203 | 42,000 | 223,000 | 0.23 x 72 | 12 | Rapid separation of lipoprotein for examination |
| | P32AT*† | 5720211204 | 32,000 | 111,000 | 12 x 32 | 186 | Simultaneous processing of many specimens |
| | P27A † | 5720211206 | 27,000 | 106,000 | 160 x 6 | 352 | Separation of large volume samples |
| P21A2 | 5720211205 | 21,000 | 71,000 | 230 x 6 | 486 | | |
| Neo Angle Rotor | P90NT*† | 5720212201 | 90,000 | 646,000 | 5 x 8 | 10 | Rapid separation of Plasmid DNA and RNA |
| | P65NT*† | 5720212202 | 65,000 | 402,000 | 12 x 10 | 23 | |
| | P65NT2*† | 5720212203 | 65,000 | 431,000 | 5 x 18 | 15 | Processing of many tubes of Plasmid DNA |
| Vertical Rotor | P100VT* | 5720213001 | 100,000 | 700,000 | 5 x 8 | 6 | Rapid separation of Plasmid DNA |
| | P65VT2*† | 5720213101 | 65,000 | 416,000 | 5 x 16 | 10 | Separation of plasmid DNA by the sedimentation equilibrium centrifugation |
| | P65VT3*† | 5720213102 | 65,000 | 402,000 | 12 x 10 | 13 | |
| | P50VT2*† | 5720213103 | 50,000 | 243,000 | 40 x 8 | 36 | Density gradient centrifugation of large volume samples |
| Swinging Bucket Rotor | P65ST † | 5720214101 | 65,000 | 419,000 | 5 x 3 | 48 | Density gradient centrifugation of Protein, DNA and RNA |
| | P56ST † | 5720214102 | 56,000 | 409,000 | 4 x 6 | 54 | Density gradient centrifugation of micro volume samples |
| | P55ST2 | 5720214001 | 55,000 | 366,000 | 5 x 6 | 50 | |
| | P40ST | 5720214002 | 40,000 | 284,000 | 13 x 6 | 139 | Density gradient centrifugation of intermediate volume samples |
| | P32ST | 5720214003 | 32,000 | 180,000 | 40 x 6 | 198 | Density gradient centrifugation of large volume samples |
| | P32ST2 | 5720214004 | 32,000 | 193,000 | 16 x 6 | 216 | Density gradient centrifugation of medium volume samples |
| | P28S | 5720214005 | 28,000 | 141,000 | 40 x 6 | 252 | Density gradient centrifugation of intermediate volume samples |
| Zonal Rotor | P32ZT | 5720215101 | 32,000 | 102,000 | 1,690ml | 363 | Large-scale purification of protein and virus particles by density gradient centrifugation |
| Continuous Flow Rotor | P32CT † | *** 5720216101 | 32,000 | 102,000 | 430ml | 42 | Continuous concentration of large virus particles |

Remarks

- * Tubes, caps and adapters are excluded in the standard accessories. Please order them separately.
- † The rotor is made-to-order item.
- ** : Rotor P50AT4 has cavities located in outer and inner. See its cross-area. K factors are 32 for outer and 38 for inner.
- Carbon fiber rotors cannot be used with CP-NX series.
- *** : To use P32CT with CP-NX or other former models, optional accessory assy is mandatory required (sold separately).
Material number of the optional accessory assy for CP-NX series is 5720410108, material number of the optional accessory assy for CP-WX and CP-MX series is 5720850525.
- Rotors with model names including the letter T are made of titanium alloy. The other rotors with names excluding the letter T are made of aluminum alloy.
All buckets for swinging bucket rotors are made of titanium alloy.
- To use zonal rotor model P32ZT, seal attachment assy model RPZ-S (Material No. 5720410109) is required. (order separately)
- When using a seal tube in the above rotors, tube sealer and respective tube rack are required. (order separately)
- Capacities of tube and bottle show their nominal. Actual capacity may not be the same to its nominal, depending on rotor structure, tube shape.





Rotor Accessories



Tube Sealer model STF3

The STF3 is heat-welding tube sealer to use the seal tubes with fixed angle rotors, vertical rotors and neo angle rotors.

Features

- ◇ Easy & Simple Usage
 1. Set seal tubes in the appropriate tube rack then set the tube rack on the rack guide of the STF3.
 2. Adjust position of the tube rack in order the tube inlet comes under a heater of the STF3.
 3. Pull the handle down to contact the heater the tube inlet and hold the handle for 1 to 2 seconds. Then pull the handle down to the end and wait until the HEAT Lamp goes off. (It takes about 40 seconds.)
- ◇ HEAT Lamp indicates the status of the heater.
- ◇ Tube racks are compatible with our former model STF2.

Specifications

| | |
|-----------------|---|
| Model | STF3 |
| Seal Method | Direct welding by the molding heater |
| Dimensions (mm) | 132 (W) x 225 (D) x 260 (H) <when raising the handle : 320 (H)> |
| Weight (kg) | 8.0 |
| Power | 110-120V version : 99 - 132 VAC, 2A, 50/60Hz 200-240V version : 180 - 264 VAC, 2A, 50/60Hz |
| Material Number | 5720310010 for 110-120V version 5720310020 for 200-240V version |

Tube Rack (sold separately)

| Material Number | Description | Tube Size |
|-----------------|----------------|-----------------|
| 5720410000 | Tube Rack (G2) | 1.5PA seal tube |
| 5720410001 | Tube Rack (G) | 2PA seal tube |
| 5720410002 | Tube Rack (B2) | 3.5PA seal tube |
| 5720410003 | Tube Rack (B3) | 4PA seal tube |
| 5720410004 | Tube Rack (B) | 5PA seal tube |
| 5720410005 | Tube Rack (B4) | 6.5PA seal tube |
| 5720410006 | Tube Rack (C2) | 8PA seal tube |
| 5720410007 | Tube Rack (C) | 12PA seal tube |
| 5720410008 | Tube Rack (E) | 40PA seal tube |
| 5720410009 | Tube Rack (F2) | 94PA seal tube |
| 5720410010 | Multi Rack* | |

* For 2, 4, 5, 6.5, 12 and 40PA seal tubes.
* Tubes sold separately as optional.


S-Cap Series

himac original S-cap series offer simple and easy operation of the tube cap for open-top (thin-walled) tubes. The unique system requires just three components and three steps to assemble the cap to the tube.

- * No consumables, such as O-ring
- * Easy cleaning

Three assembling steps

1. Fill sample into the tube up to 80% volume and insert aluminum, AL Stem into the tube by the tool.
2. Assemble titanium, Ti Ring on the tube and tighten AL Stem and Ti Ring by the tool.
3. Fill the sample into the tube by an injector and close stem by Setscrew.



Disassemble procedure after the centrifugation.

1. Remove Setscrew and withdraw 20% supernatant from the center hole of AL Stem by the injector.
2. Remove Ti Ring to downwards manually.
3. Connect tube setter into the center hole of AL Stem and remove the AL Stem by moving the tube setter left and right slowly and gradually pulling up AL Stem.

| Material Number | Descriptions | S-Cap | Tool for S-Cap | Applicable Tubes |
|-----------------|---------------------|----------|----------------|-------------------------|
| 5720410101 | S-12AL Cap Assy | 1 pair | No | 12PA tube and 12PE tube |
| 5720410102 | S-40AL Cap Assy | 1 pair | No | 40PA tube and 40PE tube |
| 5720410103 | Tool for S-Cap | No | 1 set | |
| 5720410104 | S-12AL Cap Tool Set | 8 pairs | 1 set | 12PA tube and 12PE tube |
| 5720410105 | S-12AL Cap Tool Set | 12 pairs | 1 set | 12PA tube and 12PE tube |
| 5720410106 | S-40AL Cap Tool Set | 8 pairs | 1 set | 40PA tube and 40PE tube |
| 5720410107 | S-40AL Cap Tool Set | 12 pairs | 1 set | 40PA tube and 40PE tube |

* S-cap consists of Setscrew, Al stem and Ti Ring as a pair.

* Tube excluded and sold separately.



Rotors for Density Gradient Centrifugation

The density gradient centrifugation is useful to separate multiple nano-sized particles simultaneously in the density gradient solution. We offer exclusive rotors, zonal rotor model P32ZT and continuous flow rotor model P32CT, for larger volume process.


Zonal centrifugation system : By rotor P32ZT



P32ZT
Material No. 5720215101


In zonal centrifugation, a density gradient is developed within the zonal rotor running at low speed. Next, the samples to be separated are loaded through the center of the zonal rotor. Then increase the speed to the preset high speed to separate the sample into the respective band in the density gradient in the rotor. After separation, separated sample with the density gradient solution is collected through the center by injecting the highest density gradient solution from the outside wall of the zonal rotor while the rotor is running at low speed. Collected sample with the density gradient solution is separated into the fractions. Optical density of each fraction is checked by the spectrophotometer to find the fractions containing the separated samples.

System flow diagram



- Rotor P32ZT: Excluding Seal Attachment Assy, RPZ-S, Material No. 5720410109 (sold separately).
- Density Gradient Pump required for density gradient centrifugation as optional.
- Sample Injection Pump, Density Gradient Pump and Fraction Collector not supplied from us. For detail, contact us.

Separation Procedures



Continuous flow centrifugation system : By rotor P32CT



P32CT
Material No. 5720216101


The continuous flow rotor model P32CT, designed for use with himac ultracentrifuges, enables you to perform highly efficient continuous flow separation and purification of large volume samples under high centrifugal force using differential pelleting or the density gradient centrifugation technique.

Optional 940mL Core Assy (2)

(Material No. 5720411304)


For separating samples containing much sediment or performing separation after increasing the density gradient, the optional 940mL Core Assy (2) increases the rotor capacity from 430mL to 940mL. The flow volume performance is about 35% less than that of standard core for 430mL capacity.

System flow diagram



- Rotor: P32CT rotor body with 430mL Core, Lubricating Unit, Seal Assembly, Tools (Torque wrench ,etc.), O-ring, Packing, Grease and Logbook (only P32CT with optical adapter). Seal Attachment Assy required : To use P32CT with CP-NX, the optional accessory kit is mandatory required (sold separately).
- Optional Accessory Kit exclusively for CP-NX series (Material No. 5720410108) Including: Door Assy, Rear Panel and Unit Base.
- Excluding: Sample Injection Pump, Water Cooling Bath, Density Gradient Pump and Fraction Collector. These are not supplied from us.
- Density Gradient Pump: Required for density gradient centrifugation as optional. This pump not supplied from us.

Flow rate characteristics



For details of the above, contact us. Optional Accessory Kit for CP-WX and CP-MX series, Material No. 5720850525.

Specifications

| Model | CP100NX | CP90NX | CP80NX |
|---|---|---------------------|---------------------|
| Maximum Speed | 100,000 rpm | 90,000 rpm | 80,000 rpm |
| Maximum RCF by optional rotor | 803,000 xg P100AT2 | 700,000 xg P90AT | 615,000 xg P80AT |
| Speed Control Accuracy | +/- 2rpm (1,000rpm – max. speed) | | |
| ACCEL/DECEL Mode | ACCEL : 10 / DECEL : 11 (10 and coasting) | | |
| Set Speed Range | 1,000rpm to max. speed with 100rpm increments | | |
| Timer | 1 min to 999 hours 59 min (with 1 min increments) with HOLD function | | |
| Set Temperature Range | 0°C to 40°C with 1°C increments (Accuracy : +/- 0.5°C) | | |
| Vacuum System | Oil rotary vacuum pump and oil diffusion pump | | |
| Noise Level | 51dB(A) (running at max. speed, under in-house test condition) | | |
| Heat Radiation Into a Room | 1kW or below | | |
| Cooling System | Thermo-module cooling system (CFC/HCFC/HFC-free) | | |
| Brake System | Regenerative Braking (with energy recovery) | | |
| Control Panel | Touch-sensitive color LCD panel (6.5 inch) | | |
| Data Communication | USB : Host x 1, Device x 1 / LAN x 1 | | |
| Rotor Life Management | Automatic Rotor Life Management for rotors with RLM adapter Registered by serial number to the system for rotors with optical overspeed disk | | |
| Dimensions | Width 790 x depth 690 x height 880 mm (depth 890 mm including door cover) From floor, height 925 mm to top of door handle, and 863 mm to table | | |
| Floor area | 0.81 m ² (900 x 900 mm) | | |
| Weight | 390 kg | | |
| Power Requirements | Single phase, AC200V, 208V and 220V +/-10%, max. 20A (normally 8A) Single phase, AC230V and 240V +/-10%, max. 16A (normally 7A) | | |
| Installation Environment | Ambient temperature for operation : 2°C to 40°C Ambient temperature for performance guarantee : 10°C to 30°C | | |
| Warranty of Drive Unit | 10 years after the shipment | | |
| Material No. (Asia/Pacific/Middle East) | 5720110011 | 5720109011 | 5720108011 |
| Material No. (EU market) | 5720110012 | 5720109012 | 5720108012 |
| Material No. (US/Canada) | 5720110013 | 5720109013 | 5720108013 |

Remarks :

- 1) * For detail of standard, contact us.
- 2) CP-NX series is not registered as medical device in Japan.
- 3) Due to safety reasons, installation environments, operation environments and conditions may be restricted.

- This catalogue is for international reference and not intended for a specific country.
- Orders are subject to product-availability in each country.
- All specifications are subject to change without advanced notice.
- Actual color may vary from the color of the photos on this catalogue, due to printing conditions.
- Due to safety reasons, installation environments, operating environments and conditions may be restricted.
- Unless specially mentioned, products and/or operation panel of the photos are standard specifications.
- For further information, please contact us.



CAUTION:

For safety and proper use of your machine, carefully read and follow the instruction manual.

**Eppendorf Himac Technologies Co., Ltd.
Sales Department**

2-5-12, Higashikanda, Chiyoda-ku,
Tokyo 101-0031 Japan
Tel : +81-3-5829-3616
Fax : +81-3-5829-3613
<https://www.himac-science.com/>

EHT-1C | 2022.10